- If the device plays "Pairing timeout", please try again from Step 1.
- If the device plays "Whitelist is full. Please confirm", please check the whitelist or contact after service.

To unbind EW-04 with the pendant:

Send a text message to the pendant with the format: "RBWL+ MAC address of EW-04"

Example: RBWL20:B1:7C:91:DE:B3



USER GUIDE

Check the MAC address of EW-04 on the back of the package.

Notice:

The MAC address is necessary for pairing or unbinding with the device via the configuration tool or text message. Please keep it safe.

First Connection

<u>Step 1</u>: Press and hold the side-down button + SOS button of the Pendant at the same time for 3 seconds.



Step 2: The pendant will vibrate and play a voice prompt "Start pairing, press pairing device button please, it will stop pairing after 30 seconds".



30 sec.

<u>Step 3</u>: After the voice prompt, double-click the button on the back of the EW-04.Wait the device playing a voice prompt "Pairing success".

